



Your Partner for Proteomics

ChromaTec

Recombinant Mouse Platelet Factor 4

Source	E.coli
Molecular Mass	8.2 kDa
Amino acid sequence	M – G – P – E – E etc.
Purity	> 95% as determined by SDS-PAGE and stained by Coomassie
Formulation	lyophilised from a 0.22 µm filtered solution in PBS without any detergents. Reconstitute in distilled water (our working solution for electrophoresis and N-terminal sequencing is 1 µg/µl)
Storage	lyophilised samples are storable for about twelve month Store dark at 4°C after arrival. Prevent reconstituted samples from freeze-thaw cycles. Store at -70°C after reconstitution.

www.chromatec.de

Research • Development • Production

