



Your Partner for Proteomics

ChromaTec

Recombinant Granulysin (GNLY)

Source	E. coli
Molecular Mass	18,1 kDa (15,1 kDa GNLY, double His-tagged)
Amino acid sequence	MGSSHHHHHSSGLVPRGSHMMEGL etc.
Purity	> 95% as determined by SDS-PAGE and stained by Coomassie
Formulation	lyophilised from a 0.22 µm filtered solution in A. dest with TFA without any detergens.
Storage	Store dark at 4°C after arrival. Prevent reconstituted samples from freeze-thaw cycles. Store at -70°C after reconstitution.

www.chromatec.de

Research · Development · Production

